

## Lecture 16

Cotransduction can also be used in the context of a 3-factor cross to determine gene order.

For example, say we found that two different  $Mot^-$  mutations (called  $Mot1^-$  and  $Mot2^-$ ) were found to be linked to a given Tn5 insertion each by about 50% cotransduction each (given the the error on the measurements there is no significant difference that can be used to obtain the order). To do a 3-factor cross to obtain the order of the Tn5 and  $Mot1^-$  and  $Mot2^-$  mutations, it is necessary to set up two reciprocal crosses:


### Cross 1

- (1) Grow P1 on Tn5  $Mot1^-$
- (2) Use the resulting lysate to infect a  $Mot2^-$  mutant
- (3) Select for  $Kan^r$  mutation and then count proportion that are  $Mot^+$

### Cross 2

- (1) Grow P1 on Tn5  $Mot2^-$
- (2) Use the resulting lysate to infect a  $Mot1^-$  mutant
- (3) Select for  $Kan^r$  mutation and then count proportion that are  $Mot^+$

Let's consider the crossovers that would be needed to produce  $Kan^r$   $Mot^+$  recombinants for the two possible orders.



Comparison of the frequency of  $Mot^+$  for Cross 1 and Cross 2 will immediately give us the correct order of  $Mot1^-$  and  $Mot2^-$ . For Order 1 there should be a higher frequency of  $Mot^+$  for Cross 2 than for Cross 1 (since a quadruple crossover will be required for Cross 1 but only a double crossover for Cross 2), whereas for Order 2 there should be a higher frequency of  $Mot^+$  for Cross 1 than for Cross 2.

Now let's consider how one might find a Tn5 insertion that is linked to a particular mutation (such as *MotA*<sup>-</sup>) in the first place. A simple but very laborious method would be to screen random Tn5 insertions one at a time for linkage to *MotA*<sup>-</sup> by P1 transduction. This process would require about 50 separate P1 transductions experiments. However there is a clever way to screen that would allow the desired linked Tn5 insertion to be identified in a single transduction experiment.

- 1) Start with a collection of random Tn5 insertions into wild type *E. coli* (the isolation of such a collection was described in last lecture). Grow phage P1 on the mixture of  $2 \times 10^4$  different Tn5 insertion mutants. Note that this donor strain is *MotA*<sup>+</sup>.
- 2) Use the resulting P1 phage to infect a *MotA*<sup>-</sup> recipient strain. Select for transduction of the Tn5 insertions by selecting for growth of the transductants on kanamycin plates. Screen for cotransduction of *MotA*<sup>+</sup> by testing each of the Kan<sup>r</sup> transductants for motility on soft agar. The desired cotransductant will be Kan<sup>r</sup> and will be motile. Given that one P1 phage headfull corresponds to about 1/50 of the *E. coli* chromosome, about 1 in 50 Tn5 insertions will be close enough to the *MotA* gene to show 50% cotransduction. Thus if we test about  $10^3$  Kan<sup>r</sup> transductants for motility, we are likely to find at about 10 that have cotransduced the *MotA*<sup>+</sup> marker.
- 3) Once a Tn5 (Kan<sup>r</sup>) *MotA*<sup>+</sup> transductant has been identified, P1 can be grown on this strain.
- 4) Use the P1 phage from step (3) to infect a *MotA*<sup>-</sup> recipient strain. Select for transduction of the Tn5 insertions by selecting for growth of the transductants on kanamycin plates.

Test the resulting Kan<sup>r</sup> transductants for their motility. The transductants that have cotransduced the *MotA*<sup>+</sup> marker will be motile, whereas the transductants still contain the *MotA*<sup>-</sup> allele will be nonmotile. The fraction of the total transductants that are motile will give the distance between *MotA*<sup>-</sup> and the Tn5 insertion as a cotransduction frequency.

A Tn5 (Kan<sup>r</sup>) *MotA*<sup>-</sup> transductant isolated in step (4) can then be used to transduce the *MotA*<sup>-</sup> marker into a new recipient strain by cotransduction with Tn5.